

Diamond Turbo Grinder



The diamond grinder system with the extraordinary grinding properties, extremely long service life and thus maximum efficiency due to high grade diamond and binding material. They can cut through Zirconia with a minimal amount of pressure thus alleviates the chance of micro fractures within the Zirconia substructure.





Length in mm

4.5



Order No.	TG01	TG02	TG03	TG06	TG07
Size ¢ 1/10mm	0.65	050	045	160	220
Length in mm	13.0	13.0	13.0	1.5	2.0

 Fine Grit 8000-15000 rpm max:25000 rpm 1 pc 					T
Order No.	TG08	TG09	TG10	TG11	TG12
Size ¢ 1/10mm	220	150	130	220	150

2.0


2.0

2.0

Extra Fine 8000-15000 rpm max:25000 rpm 1 pc			
Order No.	TF10	TF02	TF03
Size¢ 1/10mm	130	050	045

3.5

Diamond polishers for porcelain



3 steps diamond polishing system used to restore a natural high shine to all kind of ceramics and metals, High concentration diamond grits make more efficiency, Special formula ensure long service, odorless



fine 7000–15000 rpm max:20000 rpm 1pc		T		
Order No.	HP07F	HP06F	HP03F	HP10F
Dimensions(mm)	26*2	20*1	4*13	12*2



Separating Discs & Roughing wheels



Special dual-bonded resin binder ensure maximum durability and performance, odorless Each product has its own characteristics and is designed for specific purpose

Separating Discs (Aluminum Oxide)

"Fast-cutting!", for all non-precious metals and cobalt chrome economically, highest durability



Order No.	Dimensions	Piece
380718	38mm*0.7mm	100



Order No.DimensionsImage: Piece35071835mm*0.7mm100



Order No.	Dimensions	Piece
250618	25mm*0.6mm	100
251018	25mm*1.0mm	50

Roughing wheels (Aluminum Oxide)

Optimum processing of chrome cobalt frame via very aggressive grinding properties



Order No.	Dimensions	Piece
381018	38mm*1.0mm	50
381718	38mm*1.7mm	40



Order No.	Dimensions	Piece
351018	35mm*1.0mm	50
351718	35mm*1.7mm	40
353018	35mm*3.0mm	25

Order No.


220218




Order No.	Dimensions	Piece
322218	32mm*2.2mm	35

Separating Discs (Silicon Carbide)

For precious alloys, ceramics without leaving mark





100

22mm*0.25mm




Aluminum Oxide reinforced discs

For non-precious and partial denture alloys. Very stable even in case of lateral pressure load and high speed.





Order No.	Dimensions	Piece
400518R	40mm*0.5mm	20
400718R	40mm*0.7mm	20
401018R	40mm*1.0mm	20



Order No.	Dimensions	Piece
380718R	38mm*0.7mm	20
381018R	38mm*1.0mm	20



Order No.	Dimensions	Piece
250618R	25mm*0.6mm	20



Order No.	Dimensions	Piece
220318R	22mm*0.3mm	20

Silicon carbide reinforced discs

For separating and finishing ceramic safety




Order No.	Dimensions	Piece
220318CR	22mm*0.3mm	20



Order No.	Dimensions	Piece
220218CR	22mm*0.2mm	20

Diamond reinforced discs

Fully interspersed with diamond particles. For porcelain, Zirconia and lithium disilicate.





Order No.	Dimensions	Piece
220218DR	22mm*0.2mm	20

Diamond polishers for E.max


2 steps diamond polishing system for Lithium Disilicate from smoothing to high gloss polishing .making glazing unnecessary High concentration diamond grits make more efficiency, Special formula ensure long service, odorless





 fine 7000-15000 rpm max:20000 rpm 1pc 				
Order No.	HB07F	HB06F	HB03F	HB10F
Dimensions(mm)	26*2	20*1	4*13	12*2

Diamond polishers for zirconia







Sintered D	Diamonds					
5			-			ng
 Supercoarse 10000-1500 max:20,000r 1 pc 	0 rpm					
Order No.	SD301	SD	302 SD3	03 SI	D304	SD305
Size¢ 1/10mm	060	06	0 06	0 0	060	050
 Supercoarse 10000-1500 max:20,000r 1 pc 	0 rpm					
Order No.	SD307	SD	308 SD3	09 SI	D310	SD313
Size¢ 1/10mm	050	05	0 04	0 0	050	080
 Coarse grit 10000-15000 max:20,000rg 1 pc 					-	Ī
Order No.	SD101	SD102	SD103 SD10	4 SD105	SD107	7 SD110
Size ¢ 1/10mm	060	060	060 060	050	050	050
 Medium grit 10000-15000 max:20,000rp 1 pc 					Ī	
Order No.	SD201	SD202	SD204 SD20	5 SD207	SD210) SD213
Size ¢ 1/10mm	060	060	060 050	050	050	080
 Medium grit 10000-15000 max:20,000rg 1 pc 						
Order No.	SD222	SD	223 SD2	24 SI	D225	SD226
Size ¢ 1/10mm	016	01	6 01	6	016	016


SWIFT DI	Fine grit 10000-15000 rp max:20,000rpm	om						
	Order No.	SD001	SD002	SD003	SD004	4 SD005	5 SD006	SD007
	Size ¢ 1/10mm	060	060	060	060	050	030	050
8 9 9	10000-15000 rp max:20,000rpm	m				Ĩ		
	Order No.	SD008	SD009	SD0	10	SD011	SD012	SD013
	Size¢ 1/10mm	050	040	05)	080	040	080
*	10000-15000 rp max:20,000rpm	om						
	Order No.	SD021	SD022	SD0	23	SD024	SD025	SD026
	Size¢ 1/10mm	016	016	01	3	016	016	016
*	10000-15000 rp max:20,000rpm							
	Order No.	SD401	SD402	SD403	SD404	4 SD405	5 SD406	SD407
	Size ¢ 1/10mm	060	060	060	060	050	030	050
*	10000-15000 rp max:20,000rpm	m						
	Order No.	SD408	SD409	SD4	10	SD411	SD412	SD413
	Size¢ 1/10mm	050	040	05)	080	040	080


Flexible Diamond Discs



For the fine separation of porcelain and resin veneers. Very flexible.









Very flexible Size 1/10mm


Order No.	FD2217	FD1917
Head L in mm	0.17	0.17
coating B	Double sided full	Double sided full



Order No.	FD2215	FD1915
Head L in mm	0.15	0.15
coating B	Double sided full	Double sided full



Order No.	FD2212K	FD1912K
Head L in mm	0.12	0.12
coating B	Double sided full	Double sided full








Very flexible Size 1/10mm


Order No.	FD2217A	FD1917A
Head L in mm	0.17	0.17
coating B	Double sided 3mm	Double sided 3mm





220

Order No.	FD2215A	FD1915A
Head L in mm	0.15	0.15
coating B	Double sided 3mm	Double sided 3mm






190


Order No.	FD2212X	FD1912X
Head L in mm	0.12	0.12
coating B	Double sided full	Double sided full

Fine Grit 45um






Order No.	FD2217B	FD2222B
Head L in mm	0.17	0.22
coating B	Double	sided full






Order No.	FD1922D
Head L in mm	0.22
coating B	Double sided 3mm





FD2230A Order No. FD2222A 0.22 0.30 Head L in mm coating B Double sided 3mm






flexible Size 1/10mm


1	90
- 2	

Order No.	FD1922A	FD1930A
Head L in mm	0.22	0.30
coating B	Double sided 3mm	





Order No.	FD2222	FD2230
Head L in mm	0.22	0.30
coating B	Double sided full	






flexi	b	le	È.,		
Size	1	1	10)n	۱m

190

Order No.	FD1922	FD1930
Head L in mm	0.22	0.30
coating B	Double	sided full



Order No.	SH394	SH392
Diameter	22mm	19mm
Thickness	0.20mm	0.20mm

cleaning stone

For the cleaning and refreshing of diamond instruments particularly if you change use from metals to ceramics




Multilayer Diamond Burs



Selected natural diamond grit and advanced electroplating processing technics ensure products highest durability and sharpness











8 9	coarse 7000-15000 rpm max:25000 rpm 100pieces	-	•	
	Order No.	SP223C	SP221C	SP622C
	Dimensions(mm)	22*3	22*1.0	6*22

() () ()	medium 7000–15000 rpm max:25000 rpm 100pieces	-	-	
	Order No.	SP223M	SP221M	SP622M
	Dimensions(mm)	22*3	22*1.0	6*22

 fine 7000-15000 rpm max:25000 rpm 100pieces 	•	•	
2	0000005	000015	000005
Order No.	SP223F	SP221F	SP622F

Rubber polishers



Rubber polishers for cobalt-chrome, non-precious and semi-precious alloy



 coarse 15000-25000 rpm max:35000 rpm 100pieces 	·	\cdot	
Order No.	RM221C	RM223C	RM622C
Dimensions(mm)	22*1.0	22*3	6*22

\$ 0 0	medium 15000–25000 rpm max:35000 rpm 100pieces	•	\cdot	
	Order No.	RM221M	RM223M	RM622M
	Dimensions(mm)	22*1.0	22*3	6*22

 fine 15000-25000 rpm max:35000 rpm 100pieces 	•	•	
Order No.	RM221F	RM223F	RM622F
Dimensions(mm)	22*1.0	22*3	6*22

Silicon polishers for Acrylic



medium 5000–7000 rpm max:15000rpm 100pc				
Order No.	S1020M	S1024M	S0616M	S0517M
Dimensions(mm)	10*20	10*24	6*16	5*17


fine 5000-7000 rpm max:15000rpm 100pc				
Order No.	S1020F	S1024F	S0616F	S0517F
Dimensions(mm)	10*20	10*24	6*16	5*17

Plaster Discs



Special processing with natural diamond and stiff steel, vibration-free







Order No.	Dimensions	Piece
D853520	85mm*0.35mm*20m	ım 1



Order No.	Dimensions	9	Piece
D453518	45mm*0.35mm*1.8	mm	1

Order No.	Dimensions	9	Piece
D903020	90mm*0.30mm*20mr	m	1



Order No.	Dimensions	9	Piece
D353018	35mm*0.30mm*1.8r	nm	1

Universal polishers



For smoothing all dental materials, including porcelain, acrylics and precious metal



 coarse 7000-10000 rpm max:20000 rpm 100pieces 			
Order No.	SU515C	SU516C	SU112C
Dimensions(mm)	6*16	5*17	12*2

Mounted Point Stone



Excellent processing performance Concentric,vibration-free rotaion Long lasting



For fast processing partial denture cobalt-chrome and non-precious alloys



Dimensions (mm)	6*21	8*16	6*13
nsions (mm)	6*21	8*16	6*1

Silicon Carbide Standard grit for processing of ceramics and metal.



Made of pure aluminum oxide, or non-precious metalstal


 medium 15000-25000 rpm max:35000rpm 100pc 		T	1	Ī
Order No.	W13	w1	w20	w11
Dimensions (mm)	6*21	8*16	4.5*13	13*2.0

Made of pure aluminum oxide for precious metals and other metals

 fine 15000-25000 rp max:35000rpm 100pc 	om E		
Order No.	P45	P20	P11
Dimensions (mm)	5*17	4.5*13	13*2.0



Accessories





Mini felt wheels



 soft
 medium

 Order No.
 2504
 2505

 Dimensions (mm)
 25*4
 25*5